**Gebrauchsinformation  
User Instructions**

# **N-Bond**

**Lichthärtendes Dental-Adhäsiv  
Light curing Dental adhesive**

**CE 0124**

**MEGADENTA**

Dentalprodukte

---

MEGADENTA Dentalprodukte GmbH  
D-01454 Radeberg, Germany  
Telefon +49(0)3528- 4530, Fax +49(0)3528-45321  
[www.megadenta.de](http://www.megadenta.de) / [info@megadenta.de](mailto:info@megadenta.de)

# Gebrauchsinformation

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N-Bond

## 2. Composition:

2,2-Bis-4(2hydroxy-3-methacryloxy-propyloxy)phenyl-propane

2,2-Bis-4(2-methacryloxy-ethoxy)-phenyl-propane

3,6-Dioxaoctamethylendimethacrylat

Hydroxypropylmethacrylat

Nano-filler , Initiators

## 3. Indication

Dental adhesive for restorative dentistry, repairs of acrylic-, ceramic ore veneering parts.

## 4. Contraindications

So far unknown

## 5. Interactions

N-Bond may not be used together with under fillings containing Eugenol because Phenolic substances could inhibit the polymerization process.

## 6. Application method:

- **Premedication:** Prepare the cavity and dry it (rubberdam), do not remove the smear layer and clean unprepared dentine with a cleaning paste. The prepared cavity is now rinsed with water and dried with oil free air. Deeper lying areas of the cavity can be filled with a suitable under filling material.

- **Enamel and dentine etching:** The cleaned enamel edges are etched with C-Cid. C-Cid is applied with the supplied brush. After 20-30 s C-Cid is applied on the dentine allowed to take effect for another 10-15 s and then thoroughly rinsed with water. Excessive drying should be avoided in this case as this would cause the exposed collagen fiber to collapse.

-**Application:** One drop of N-Bond is dispensed onto the small brush and applied as a thin coating to the etched surface and spreaded evenly with a gentle flow of air. Now the composite can be applied immediately or by choice, after curing (10s) with a blue light source (e.a. Megalux CS) of the bonding layer. A contamination of the bonded surface with blood or saliva has to be avoided. If not, the surface has to be re-etched again.

Following the polymerization of N-Bond, a light curing composite (e.g. N-Fill / Megafill MH) can be applied. Then proceed as instructed by the composite manufacturer.

## 7. Storage stability

3 years

## 8. Storage and keeping instructions

Refrigerated storage (+5°C to +8°C) is recommended; do not exceed 25°C.

## 9. Administration and package size

5,0 ml N-Bond bottle as single package ore part of assortment packages

## 10. Date information published

07/2008

## 11. Name of manufacture

MEGADENTA Dentalprodukte GmbH, D-01454 Radeberg, Germany

Tel. +49-3528-453-0 Fax +49-3528-453-21 Mail : info@megadenta.de